

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 101

**Grad. pipettes, 1:0,01 ml, AS, conf. cert.**




**Batch:** 07-19

**Serial:**

**Measuring Nr.** 28098

**Batch size:** 4296

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	0,9982 ml	0,5985 ml	0,2002 ml
Deviation:	0,0022 ml	0,0021 ml	0,0012 ml
Precision (CV):	0,2156 %	0,3487 %	0,5898 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	04.03.2019
Balance:	Sartorius - LE244S-OCE

