

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 101

Grad. pipettes, 1:0,01 ml, AS, conf. cert.




Batch: 15-20

Serial:

Measuring Nr. 37376

Batch size: 2459

| Volume: | 1,0000 ml | 0,6000 ml | 0,2000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0070 ml | ± 0,0070 ml | ± 0,0070 ml |
| Average value: | 0,9990 ml | 0,6016 ml | 0,2008 ml |
| Deviation: | 0,0020 ml | 0,0022 ml | 0,0013 ml |
| Precision (CV): | 0,2005 % | 0,3721 % | 0,6589 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 9 |
| Date: | 28.04.2020 |
| Balance: | Sartorius - LE244S-OCE |

