

# Certificate of performance



**Article:** 3 510 110 BL

**Grad. pipettes, 1:0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 08-18

**Serial:**

**Measuring Nr.** 19066

**Batch size:** 280

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
<b>Tolerance:</b>	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
<b>Average value:</b>	0,9980 ml	0,5970 ml	0,1989 ml
<b>Deviation:</b>	0,0020 ml	0,0010 ml	0,0005 ml
<b>Precision (CV):</b>	0,1999 %	0,1692 %	0,2301 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	44
<b>Date:</b>	23.03.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

