

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 110

**Grad. pipettes, 1:0,1 ml, AS, conf. cert.**

**Batch:** 09-23




**Serial:**

**Measuring Nr.** 61871

**Batch size:**

280

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	0,9986 ml	0,5997 ml	0,1999 ml
Deviation:	0,0023 ml	0,0007 ml	0,0008 ml
Precision (CV):	0,2300 %	0,1158 %	0,4002 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

02.03.2023

**Balance:**

Sartorius -  
LE244S-OCE

DE-M