

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 202 BL

**Grad. pipettes, 2:0,02 ml, blue grad., AS, conf. cert.**




**Batch:** 05-24

**Serial:**

**Measuring Nr.** 65505

**Batch size:** 280

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9991 ml	1,2035 ml	0,4021 ml
Deviation:	0,0022 ml	0,0024 ml	0,0011 ml
Precision (CV):	0,1088 %	0,2016 %	0,2790 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	06.02.2024
Balance:	Sartorius - LE244S-OCE

DE-M