

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 202 USP

**Grad. pipettes, 2:0, 0,2 ml, USP, AS, conf. cert.**




**Batch:** 40-23

**Serial:**

**Measuring Nr.** 64781

**Batch size:** 150

| Volume:         | 2,0000 ml   | 1,2000 ml   | 0,4000 ml   |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value:  | 1,9996 ml   | 1,1993 ml   | 0,4007 ml   |
| Deviation:      | 0,0023 ml   | 0,0024 ml   | 0,0011 ml   |
| Precision (CV): | 0,1162 %    | 0,2023 %    | 0,2800 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                           |
|-----------|---------------------------|
| Operator: | 7                         |
| Date:     | 24.10.2023                |
| Balance:  | Sartorius -<br>LE244S-OCE |

DE-M