

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 202

Grad. pipettes, 2:0,02 ml, AS, conf. cert.




Batch: 41-17

Serial:

Measuring Nr. 15183

Batch size: 1029

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	2,0005 ml	1,1980 ml	0,4002 ml
Deviation:	0,0034 ml	0,0024 ml	0,0014 ml
Precision (CV):	0,1698 %	0,2031 %	0,3587 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	19.10.2017
Balance:	Sartorius - LE244S-OCE

