

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 202

Grad. pipettes, 2:0,02 ml, AS, conf. cert.




Batch: 25-22

Serial:

Measuring Nr. 56755

Batch size: 2983

| Volume: | 2,0000 ml | 1,2000 ml | 0,4000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value: | 2,0010 ml | 1,2006 ml | 0,4019 ml |
| Deviation: | 0,0036 ml | 0,0029 ml | 0,0016 ml |
| Precision (CV): | 0,1813 % | 0,2377 % | 0,4018 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 7 |
| Date: | 24.06.2022 |
| Balance: | Sartorius - LE244S-OCE |

