

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 202

**Grad. pipettes, 2:0,02 ml, AS, conf. cert.**

**Batch:** 25-23




**Serial:**

**Measuring Nr.** 63541

**Batch size:**

3118

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9993 ml	1,2024 ml	0,4014 ml
Deviation:	0,0023 ml	0,0024 ml	0,0014 ml
Precision (CV):	0,1171 %	0,1991 %	0,3478 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

28.06.2023

**Balance:**

Sartorius -  
LE244S-OCE

