

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 202

Grad. pipettes, 2:0,02 ml, AS, conf. cert.




Batch: 43-23

Serial:

Measuring Nr. 64894

Batch size: 3200

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9995 ml	1,1989 ml	0,4006 ml
Deviation:	0,0034 ml	0,0025 ml	0,0016 ml
Precision (CV):	0,1684 %	0,2108 %	0,3923 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	31.10.2023
Balance:	Sartorius - LE244S-OCE

