

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 202

**Grad. pipettes, 2:0,02 ml, AS, conf. cert.**




**Batch:** 14-24

**Serial:**

**Measuring Nr.** 66735

**Batch size:** 3199

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	1,9993 ml	1,2014 ml	0,4012 ml
Deviation:	0,0033 ml	0,0028 ml	0,0018 ml
Precision (CV):	0,1636 %	0,2344 %	0,4410 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	11.04.2024
Balance:	Sartorius - LE244S-OCE

