

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 210

**Grad. pipettes, 2:0,1 ml, AS, conf. cert.**




**Batch:** 39-17

**Serial:** -

**Measuring Nr.** 14809

**Batch size:** 500

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
<b>Tolerance:</b>	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
<b>Average value:</b>	1,9981 ml	1,1975 ml	0,3999 ml
<b>Deviation:</b>	0,0027 ml	0,0018 ml	0,0018 ml
<b>Precision (CV):</b>	0,1368 %	0,1515 %	0,4524 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

11.10.2017

**Balance:**

Sartorius -  
LE244S-OCE

DE-M