

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 505 BL USP

**Grad. pipettes, 5:0,05 ml, blue grad., AS, conf.cert.**




**Batch:** 05-19

**Serial:**

**Measuring Nr.** 27513

**Batch size:** 150

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
<b>Average value:</b>	4,9940 ml	2,9984 ml	0,9982 ml
<b>Deviation:</b>	0,0056 ml	0,0038 ml	0,0025 ml
<b>Precision (CV):</b>	0,1130 %	0,1260 %	0,2462 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	06.02.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

