

Certificate of performance



Article: 3 510 505 USP

graduated pipettes 5 ml class AS USP KB

Batch: 08-22




Serial:

Measuring Nr.: 52151

Batch size:

280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	4,9914 ml	3,0000 ml	1,0003 ml
Deviation:	0,0052 ml	0,0035 ml	0,0023 ml
Precision (CV):	0,1030 %	0,1173 %	0,2322 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	17.02.2022
Balance:	Sartorius - LE244S-OCE

