

Certificate of performance



Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 36-20

Serial:

Measuring Nr.: 40067

Batch size: 3130

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9894 ml	2,9950 ml	0,9993 ml
Deviation:	0,0043 ml	0,0043 ml	0,0037 ml
Precision (CV):	0,0859 %	0,1432 %	0,3705 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	15.09.2020
Balance:	Sartorius - LE244S-OCE

