

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.

Batch: 19-22




Serial:

Measuring Nr. 55542

Batch size:

3032

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9977 ml	2,9991 ml	1,0029 ml
Deviation:	0,0039 ml	0,0046 ml	0,0031 ml
Precision (CV):	0,0782 %	0,1537 %	0,3106 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

19.05.2022

Balance:

Sartorius -
LE244S-OCE

DE-M