

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 41-22

Serial:

Measuring Nr. 59759

Batch size: 3200

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	5,0012 ml	2,9960 ml	0,9984 ml
Deviation:	0,0039 ml	0,0042 ml	0,0018 ml
Precision (CV):	0,0787 %	0,1385 %	0,1766 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	21.10.2022
Balance:	Sartorius - MSE224S-0CE -DU

