

Certificate of performance



Article: 3 510 505

Graduated pipettes, 5:0,05 ml, AS, conf. cert.




Batch: 32-23

Serial:

Measuring Nr. 64156

Batch size: 4996

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	4,9970 ml	2,9951 ml	1,0007 ml
Deviation:	0,0050 ml	0,0045 ml	0,0019 ml
Precision (CV):	0,1004 %	0,1484 %	0,1912 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	29.08.2023
Balance:	Sartorius - MSE224S-0CE -DU

