

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 510 BL

**Grad. pipettes, 5:0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 37-20

**Serial:**

**Measuring Nr.** 40066

**Batch size:** 417

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	4,9942 ml	2,9976 ml	1,0005 ml
<b>Deviation:</b>	0,0040 ml	0,0039 ml	0,0027 ml
<b>Precision (CV):</b>	0,0795 %	0,1312 %	0,2699 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	11.09.2020
<b>Balance:</b>	Sartorius - LE244S-OCE

