

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 510

Grad. pipettes, 5:0,1 ml, AS, conf. cert.




Batch: 40-22

Serial: -

Measuring Nr. 59694

Batch size: 1200

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	5,0010 ml	2,9958 ml	0,9969 ml
Deviation:	0,0042 ml	0,0036 ml	0,0017 ml
Precision (CV):	0,0842 %	0,1205 %	0,1673 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.10.2022

Balance:

Sartorius -
MSE224S-0CE
-DU

