

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710 BL USP

**Grad.pipette, 10: 0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 05-19

**Serial:**

**Measuring Nr.** 27514

**Batch size:** 150

Volume:	10,000 ml	6,000 ml	2,000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	9,9944 ml	6,0052 ml	2,0061 ml
<b>Deviation:</b>	0,0051 ml	0,0079 ml	0,0039 ml
<b>Precision (CV):</b>	0,0508 %	0,1323 %	0,1972 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	06.02.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

