

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710 USP

**Graduated pipettes, 10:0,1 ml, USP, AS, conf.cert.**




**Batch:** 25-23

**Serial:**

**Measuring Nr.** 63600

**Batch size:** 280

Volume:	10,000 ml	6,0000 ml	2,0000 ml
<b>Tolerance:</b>	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
<b>Average value:</b>	9,9990 ml	5,9892 ml	1,9992 ml
<b>Deviation:</b>	0,0082 ml	0,0046 ml	0,0033 ml
<b>Precision (CV):</b>	0,0815 %	0,0762 %	0,1660 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	06.07.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

