

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710

**Graduated pipettes, 10:0,1 ml, AS, conf. cert.**




**Batch:** 04-18

**Serial:**

**Measuring Nr.** 17916

**Batch size:** 5000

Volume:	10,000 ml	6,0000 ml	2,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	9,9892 ml	5,9920 ml	2,0014 ml
<b>Deviation:</b>	0,0085 ml	0,0090 ml	0,0044 ml
<b>Precision (CV):</b>	0,0851 %	0,1493 %	0,2214 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	13.02.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

