

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 15-18

Serial:

Measuring Nr. 19974

Batch size: 5000

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9953 ml	5,9956 ml	2,0019 ml
Deviation:	0,0074 ml	0,0115 ml	0,0042 ml
Precision (CV):	0,0742 %	0,1915 %	0,2085 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	18.04.2018
Balance:	Sartorius - LE244S-OCE

