

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710

**Graduated pipettes, 10:0,1 ml, AS, conf. cert.**




**Batch:** 27-18

**Serial:**

**Measuring Nr.** 22070

**Batch size:** 4927

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	10,0003 ml	6,0134 ml	2,0062 ml
Deviation:	0,0069 ml	0,0081 ml	0,0036 ml
Precision (CV):	0,0689 %	0,1358 %	0,1794 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	44
Date:	10.07.2018
Balance:	Sartorius - LE244S-OCE

