

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 32-19

Serial:

Measuring Nr. 32129

Batch size: 5000

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9868 ml	5,9934 ml	1,9995 ml
Deviation:	0,0071 ml	0,0085 ml	0,0046 ml
Precision (CV):	0,0710 %	0,1424 %	0,2307 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	16.08.2019
Balance:	Sartorius - LE244S-OCE

