

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.

Batch: 42-19




Serial:

Measuring Nr. 33870

Batch size:

8000

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9893 ml	5,9923 ml	2,0008 ml
Deviation:	0,0093 ml	0,0093 ml	0,0056 ml
Precision (CV):	0,0929 %	0,1553 %	0,2785 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

29.10.2019

Balance:

Sartorius -
MSE224S-0CE
-DU

