

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 14-20

Serial:

Measuring Nr. 37262

Batch size: 3200

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9865 ml	5,9904 ml	2,0057 ml
Deviation:	0,0080 ml	0,0088 ml	0,0061 ml
Precision (CV):	0,0802 %	0,1472 %	0,3073 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.04.2020
Balance:	Sartorius - MSE224S-0CE -DU

