

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 24-20

Serial:

Measuring Nr. 38813

Batch size: 5000

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9954 ml	5,9923 ml	2,0053 ml
Deviation:	0,0087 ml	0,0106 ml	0,0049 ml
Precision (CV):	0,0872 %	0,1760 %	0,2446 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	24.06.2020
Balance:	Sartorius - LE244S-OCE

