

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 35-22

Serial:

Measuring Nr. 59578

Batch size: 4999

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9929 ml	6,0000 ml	1,9991 ml
Deviation:	0,0078 ml	0,0090 ml	0,0050 ml
Precision (CV):	0,0782 %	0,1504 %	0,2478 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

06.10.2022

Balance:

Sartorius -
LE244S-OCE

DE-M