

Certificate of performance



Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.

Batch: 25-23




Serial:

Measuring Nr. 63604

Batch size:

3200

Volume:	10,000 ml	6,0000 ml	2,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9942 ml	6,0023 ml	2,0011 ml
Deviation:	0,0104 ml	0,0085 ml	0,0062 ml
Precision (CV):	0,1038 %	0,1417 %	0,3077 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	10.07.2023
Balance:	Sartorius - LE244S-OCE

