

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 710

**Graduated pipettes, 10:0,1 ml, AS, conf. cert.**




**Batch:** 30-23

**Serial:**

**Measuring Nr.** 63885

**Batch size:** 8999

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
<b>Tolerance:</b>	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
<b>Average value:</b>	9,9988 ml	6,0043 ml	2,0025 ml
<b>Deviation:</b>	0,0071 ml	0,0066 ml	0,0045 ml
<b>Precision (CV):</b>	0,0712 %	0,1108 %	0,2260 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	02.08.2023
<b>Balance:</b>	Sartorius - LE244S-OCE

