

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 710

Graduated pipettes, 10:0,1 ml, AS, conf. cert.




Batch: 26-24

Serial:

Measuring Nr. 67962

Batch size: 3200

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	9,9934 ml	5,9897 ml	1,9952 ml
Deviation:	0,0080 ml	0,0094 ml	0,0060 ml
Precision (CV):	0,0795 %	0,1572 %	0,3016 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	16.07.2024
Balance:	Sartorius - LE244S-OCE

