

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 810 BL

Grad. pipettes, 20:01, blue grad., AS, conf. cert




Batch: 44-19

Serial:

Measuring Nr. 34047

Batch size: 280

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	19,9903 ml	11,9940 ml	4,0010 ml
Deviation:	0,0097 ml	0,0184 ml	0,0100 ml
Precision (CV):	0,0486 %	0,1533 %	0,2509 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	04.11.2019
Balance:	Sartorius - LE244S-OCE

