

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 810 BL

**Grad. pipettes, 20:01, blue grad., AS, conf. cert**




**Batch:** 14-22

**Serial:**

**Measuring Nr.** 54371

**Batch size:** 280

Volume:	20,000 ml	12,000 ml	4,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	20,0133 ml	12,0008 ml	4,0112 ml
<b>Deviation:</b>	0,0149 ml	0,0175 ml	0,0075 ml
<b>Precision (CV):</b>	0,0744 %	0,1461 %	0,1870 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	08.04.2022
<b>Balance:</b>	Sartorius - LE244S-OCE

