

Certificate of performance



Article: 3 510 810

Graduated pipettes, 20 ml, AS, conf. cert.




Batch: 26-17

Serial: -

Measuring Nr.: 11759

Batch size: 1200

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	19,9796 ml	11,9970 ml	4,0015 ml
Deviation:	0,0117 ml	0,0142 ml	0,0060 ml
Precision (CV):	0,0585 %	0,1186 %	0,1489 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

44

Date:

04.07.2017

Balance:

Sartorius -
LE244S-OCE

