

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 810

Graduated pipettes, 20:0,1 ml, AS, conf. cert.




Batch: 46-18

Serial:

Measuring Nr. 25939

Batch size: 1200

| Volume: | 20,000 ml | 12,000 ml | 4,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,1000 ml | ± 0,1000 ml | ± 0,1000 ml |
| Average value: | 19,9952 ml | 12,0104 ml | 4,0109 ml |
| Deviation: | 0,0140 ml | 0,0137 ml | 0,0072 ml |
| Precision (CV): | 0,0699 % | 0,1143 % | 0,1788 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 29.11.2018 |
| Balance: | Sartorius - LE244S-OCE |

