

Certificate of performance



Article: 3 510 810

Graduated pipettes, 20:0,1 ml, AS, conf. cert.




Batch: 19-19

Serial:

Measuring Nr. 29702

Batch size: 1200

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	19,9955 ml	12,0060 ml	4,0080 ml
Deviation:	0,0142 ml	0,0168 ml	0,0082 ml
Precision (CV):	0,0709 %	0,1403 %	0,2049 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	08.05.2019
Balance:	Sartorius - LE244S-OCE

