

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 810

Graduated pipettes, 20:0,1 ml, AS, conf. cert.




Batch: 44-20

Serial:

Measuring Nr. 41835

Batch size: 1200

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	19,9833 ml	11,9851 ml	4,0067 ml
Deviation:	0,0142 ml	0,0148 ml	0,0115 ml
Precision (CV):	0,0711 %	0,1237 %	0,2878 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

28.10.2020

Balance:

Sartorius -
LE244S-OCE

DE-M