

# Certificate of performance



**Article:** 3 510 810

**Graduated pipettes, 20:0,1 ml, AS, conf. cert.**




**Batch:** 05-23

**Serial:**

**Measuring Nr.:** 61702

**Batch size:** 1200

Volume:	20,000 ml	12,000 ml	4,000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	20,0125 ml	12,0185 ml	4,0098 ml
Deviation:	0,0107 ml	0,0125 ml	0,0079 ml
Precision (CV):	0,0534 %	0,1039 %	0,1966 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	15.02.2023
Balance:	Sartorius - LE244S-OCE

