

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 910 BL

**Grad. pipettes, 25:0,1 ml, blue grad., AS, conf. cert.**




**Batch:** 49-19

**Serial:**

**Measuring Nr.** 34941

**Batch size:** 280

Volume:	25,000 ml	15,000 ml	5,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	24,9979 ml	15,0100 ml	5,0148 ml
<b>Deviation:</b>	0,0154 ml	0,0183 ml	0,0104 ml
<b>Precision (CV):</b>	0,0614 %	0,1221 %	0,2080 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	05.12.2019
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

