

Certificate of performance



Article: 3 510 910

Grad. pipettes, 25:0,1 ml, AS, conf. cert.




Batch: 40-17

Serial:

Measuring Nr. 15237

Batch size: 1200

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	24,9924 ml	14,9954 ml	5,0045 ml
Deviation:	0,0185 ml	0,0142 ml	0,0104 ml
Precision (CV):	0,0741 %	0,0945 %	0,2087 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	20.10.2017
Balance:	Sartorius - LE244S-OCE

