

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 910

**Grad. pipettes, 25:0,1 ml, AS, conf. cert.**




**Batch:** 27-18

**Serial:**

**Measuring Nr.** 22230

**Batch size:** 1200

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	24,9932 ml	14,9918 ml	4,9976 ml
<b>Deviation:</b>	0,0157 ml	0,0162 ml	0,0080 ml
<b>Precision (CV):</b>	0,0626 %	0,1082 %	0,1595 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 44

**Date:** 13.07.2018

**Balance:** Sartorius -  
LE244S-OCE

DE-M