

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 910

**Grad. pipettes, 25:0,1 ml, AS, conf. cert.**




**Batch:** 07-19

**Serial:**

**Measuring Nr.** 27659

**Batch size:** 1200

Volume:	25,000 ml	15,000 ml	5,000 ml
<b>Tolerance:</b>	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
<b>Average value:</b>	24,9950 ml	14,9864 ml	4,9999 ml
<b>Deviation:</b>	0,0133 ml	0,0187 ml	0,0090 ml
<b>Precision (CV):</b>	0,0533 %	0,1247 %	0,1792 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	14.02.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

