

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 910

Grad. pipettes, 25:0,1 ml, AS, conf. cert.




Batch: 40-19

Serial:

Measuring Nr. 33574

Batch size: 1458

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	24,9876 ml	15,0036 ml	5,0036 ml
Deviation:	0,0160 ml	0,0164 ml	0,0086 ml
Precision (CV):	0,0640 %	0,1096 %	0,1723 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	10.10.2019
Balance:	Sartorius - LE244S-OCE

