

# Certificate of performance



**Article:** 3 510 910

**Grad. pipettes, 25:0,1 ml, AS, conf. cert.**




**Batch:** 25-24

**Serial:**

**Measuring Nr.** 67793

**Batch size:** 500

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	24,9762 ml	15,0004 ml	5,0007 ml
Deviation:	0,0126 ml	0,0195 ml	0,0137 ml
Precision (CV):	0,0504 %	0,1299 %	0,2747 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	27.06.2024
Balance:	Sartorius - LE244S-OCE

