

Certificate of performance



Article: 3 510 910

Grad. pipettes, 25:0,1 ml, AS, conf. cert.




Batch: 04-25

Serial:

Measuring Nr. 71812

Batch size: 365

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	24,9860 ml	14,9861 ml	5,0000 ml
Deviation:	0,0169 ml	0,0155 ml	0,0103 ml
Precision (CV):	0,0676 %	0,1036 %	0,2053 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	29.01.2025
Balance:	Sartorius - LE244S-OCE

