

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 101

Grad. pipette, 1 : 0,01 ml, AS, conf. cert.

Batch: 35-21




Serial:

Measuring Nr. 48626

Batch size:

500

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	$\pm 0,0070$ ml	$\pm 0,0070$ ml	$\pm 0,0070$ ml
Average value:	0,9996 ml	0,6006 ml	0,2005 ml
Deviation:	0,0017 ml	0,0012 ml	0,0008 ml
Precision (CV):	0,1720 %	0,2081 %	0,3803 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.09.2021

Balance:

Sartorius -
LE244S-OCE

DE-M