

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 101

Grad. pipette, 1 : 0,01 ml, AS, conf. cert.




Batch: 31-23

Serial:

Measuring Nr. 63893

Batch size: 500

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	0,9984 ml	0,6003 ml	0,2004 ml
Deviation:	0,0022 ml	0,0017 ml	0,0007 ml
Precision (CV):	0,2203 %	0,2780 %	0,3612 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	09.08.2023
Balance:	Sartorius - LE244S-OCE

