

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 110

**Grad. pipette, 1:0,1 ml, AS, conf. cert.**




**Batch:** 49-20

**Serial:**

**Measuring Nr.** 43108

**Batch size:** 150

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	1,0011 ml	0,5986 ml	0,1998 ml
Deviation:	0,0023 ml	0,0018 ml	0,0010 ml
Precision (CV):	0,2277 %	0,3081 %	0,5073 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	77
Date:	08.12.2020
Balance:	Sartorius - MSE224S-0CE -DU

